# EYE Care Solution WIBE B1

(Color: Silver, Pink)

## Description:

EYE Care Solution WIBE B1 is a cutting-edge technology device where a blended effect of Heat and special light waves are deployed to regain the health and vitality of the under-eyes area overcoming aging marks at this sensitive area.

A heated eye mask activates the blood circulation and expands the oily eyelids glands to facilitate the discharge of wastes.

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## Max 38°C Heating procedure

Delivers a soothing warmth temperature of 38°C close to the body temperature, relaxing the muscle around the eyes and alleviating fatigue

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## LED lightening

By using 48 near Infrared rays (830 nm) and visible light (630 nm) simultaneously, Collagen activity and blood circulation are stimulated

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## Easy operation with one bottom

Wireless operation makes it easy and convenient to use anytime, anywhere, and supports 3 modes with one button

Heating & Lighting

Heating

Lighting

## How to use

1. Press and hold the power button to turn it on
2. press the button to select the desired mode.
3. Select a mode and put the mask ON by attaching it lightly to the face.
4. When taken off, the sensor stops the operation automatically.
5. resumes operation when worn again.
6. Shuts down after 15 minutes of use.

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